

# Wireless Foot Control

## 360° of Freedom

The Wireless Foot Control from ArthroCare employs invisible infrared technology to eliminate the constraints and hazards of cables associated with traditional foot controls. Because it uses light instead of radio waves, the wireless foot control resists interference with other equipment in the operating room for consistent, reliable operation.



### Wireless Foot Control

Designed to work with the Quantum™ System only.

- Easy to use
- No calibration required
- 360° Signal Coverage
- Covers up to 40 feet of signal range
- Over 160 hours battery life
- Uses AA standard batteries

## Benefits

- Easy to use, no calibration required – simply plug the receiver into the Quantum controller and look for the illuminated green indicator light
- No cumbersome cords that could cause stumbling or restrict user movement
- Infrared light is consistent and reliably covers up to 40 feet with 360 degrees of signal
- Over 160 hours of battery life under standard operating conditions\* ensures dependable performance time after time
- Uses 3 standard AA batteries
- No need to worry about interrupting a case or wasting valuable OR time; Low Battery indicator allows ample time for battery replacement after illuminating\*
- No hassle wipe down cleaning
- Designed to meet UL60601.1. CE, and EMC requirements



Quantum Controller



Wireless Foot Control and Receiver

## Technical Specifications

| Catalog Number:       | H4001-01      |                     |              |
|-----------------------|---------------|---------------------|--------------|
|                       | Size          | Weight              | Battery Life |
| Wireless Foot Control | 14"x10"x3"    | 5.5 lbs w/batteries | >160 hours*  |
| Receiver              | 5"x3.5"x2.75" | .45 lbs             | N/A          |



**ArthroCare Sports Medicine**  
 680 Vaqueros Avenue  
 Sunnyvale, CA 94085-3523  
 order entry **800-797-6520**  
 phone 408-736-0224  
 order entry fax 888-994-2782  
[arthrocarsportsmedicine.com](http://arthrocarsportsmedicine.com)

**ArthroCare Europe AB**  
 Baggensgatan 25  
 111 31 Stockholm  
 Sweden  
 phone +46 8 546 172 00  
 fax +46 8 546 172 39  
[info@arthrocare.se](mailto:info@arthrocare.se)

\*Data on file at ArthroCare Corporation, Sunnyvale, CA, USA.

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.  
 ©2007 ArthroCare Corporation. ArthroCare, Quantum and Multi-Electrode Technology are trademarks and/or registered trademarks of ArthroCare Corporation. ArthroWand surgical devices are sold exclusively by ArthroCare Corporation.

P/N 30000 Rev A